Improving Global TMJ Health, Daily!®

TMJ Metal-on-Metal Total Joint Replacement Prosthesis Systems™
comprised of the TMJ Fossa Eminence and Condylar Prostheses

Instructions for Use

TMJ Fossa Eminence Prosthesis Systems™
Right Side System
(Models FER-01 through FER-44)

Left Side System
(Models FEL-01 through FEL-44)

TMJ Patient-Specific™ Fossa-Eminence Prosthesis System
(Models CFER and CFEL)

TMJ Condylar Prostheses
TMJ Universal ArthroChrome™ Condylar Prosthesis System
(Models R/L-MCP-45, R/L-MCP-50, R/L-MCP-55)

TMJ Patient-Specific™ Condylar Prosthesis System
(Models CRMCP and CLMCP)

TMJ Christensen/Chase ArthroChrome™ Condylar Prosthesis System
(Models C/C-RMCP-45, C/C-RMCP-50, C/C-RMCP-55,
C/C-LMCP-45, C/C-LMCP-50, C/C-LMCP-55)

Store between 10° and 32°C (50° - 90°F) and 20% to 80% relative humidity.

Manufactured by
TMJ Implants, Inc.
17301 West Colfax Ave, #135
Golden, Colorado 80401
USA
303- 277-1338
800-825-4TMJ (800-825-4865)
303-277-1421 (facsimile)
email: info@tmj.com

European Authorized Representative ("E.A.R."):
Obelis s.a. (O.E.A.R.C.)
Avenue de Tervuren, 34 b44
1040 Brussels
Tel. : +32.2.732.5954
Fax : +32.2.732.6003
Email mail@obelis.net

Caution: United States Federal Law restricts this device to sale by or on the order of a physician.

1. DEVICE DESCRIPTION

The TMJ Metal-on-Metal Total Joint Replacement System™ is comprised of the TMJ Fossa-Eminence Prosthesis System™ and the TMJ Condylar Prosthesis System™, respectively.

THE TMJ FOSSA-EMINENCE PROSTHESIS

The TMJ Fossa-Eminence Prosthesis is designed to provide a thin, rigid, well-fitting prosthetic covering for the articulating surface of the temporomandibular joint comprised of the glenoid fossa and articular eminence of the temporal bone. The articular surface of the prosthetic glenoid fossa and articular eminence is highly polished to minimize friction in joint movement.

The prosthesis, and the screws with which it is to be secured to the skull, are manufactured from surgical Co-Cr-Mo alloy (ASTM F75/ASTM F799). These devices are intended for permanent implant and are for single use only.

All components in this Fossa System, including individual prosthesis, drills and screws are sterilized by gamma-irradiation or e-beam irradiation (25 kGy), and are packaged in individual double-peel PETG and Tyvek containers.

Additional **NON-STERILE** Fossa Kits or Systems accompany this TMJ Fossa-Eminence Prosthesis System™ and are essential for its use. The Fossa-Eminence Trial Sizing System contains trial sizer components for each size of implant. The Instrument Kit contains screwdrivers and Fossa-Eminence holders. (The Sizing System and Instrument Kit may be packaged together.)

The accompanying kits or systems must be steam sterilized prior to use in accordance with procedures outlined in Section 9.

THE TMJ CONDYLAR PROSTHESIS

The TMJ Condylar prosthesis, which is intended to be used with the Fossa Eminence for total joint replacement, is designed to replace the articular surface of the mandibular condyle.

The TMJ Condylar Prosthesis Systems are designed to seat against the TMJ Fossa-Eminence prosthesis and to be secured to the ramus of the mandible with cobalt chrome alloy screws. The Universal Condylar prosthesis is manufactured in three lengths, and is designed to be used on either the right or left side. The Christensen/Chase Condylar prosthesis is also available in three lengths and is manufactured specifically for either the right or left side.

The entire ArthroChrome™ prosthesis and the screws, with which they are to be secured to the mandible, are manufactured from surgical Co-Cr-Mo alloy (ASTM F75/ASTM F799). These devices are intended for permanent implant and are for single use only.

All components in the Condylar Prosthesis System, including individual prosthesis, drills and screws are sterilized by gamma-irradiation or e-beam irradiation (2.5 Mrads), and are packaged in individual double-peel PETG and Tyvek containers.

Additional Condylar kits or systems accompany this TMJ Condylar Prosthesis System and are essential for its use. The Condylar Trial Sizing System contains sterile disposable trial sizer components for each size of implant, as well as a sizing template. The **NON-STERILE** Instrument Kit or System contains screwdrivers and holders. (The Instrument Kit may be packaged with the Fossa-Eminence Sizing System.)

The accompanying instrument kit or system must be steam sterilized prior to use in accordance with procedures outlined in Section 9.

2. INTENDED USE/INDICATIONS

The TMJ Metal-on-Metal Total Joint Replacement System™ is intended for use in treatment of severe temporomandibular joint disease. The TMJ Condylar prosthesis is intended for use in conjunction with the TMJ Fossa-Eminence prosthesis whenever total joint reconstruction is necessary due to:

- Inflammatory arthritis involving the temporomandibular joint not responsive to other modalities of treatment,
- Recurrent fibrous and or bony ankylosis not responsive to other modalities of treatment,
- Failed tissue graft,
- Failed alloplastic joint reconstruction,
- Loss of vertical mandibular height and/or occlusal relationship due to bone resorption, trauma, developmental abnormality, or pathologic lesion.

3. CONTRAINDICATIONS

The TMJ Metal-on-Metal Total Joint Replacement System™ should not be used for patients with one or more of the following conditions:

- Infection or malignancy in the head or neck region,
- Known allergy to any of the components of the system,
- Ability to exert significant post-operative masticatory muscle hyperfunction (clenching or grinding), which may lead to overload and fracture of the device or loosening of the screws.

4. WARNINGS

- **READ ALL ACCOMPANYING LABELING**

Prior to use, the surgeon must read the entire Instructions for Use and device labeling.

Dynamic fatigue tests were conducted on the TMJ Implants Metal-on-Metal Total Joint Replacement System with a force applied vertically to the device. No failures occurred at 130 lbs. Physicians should carefully consider the results of these fatigue tests when patients present with particular anatomical considerations or unusual masticatory forces.

- **PERFECT THE TECHNIQUE FOR IMPLANTATION**

It is strongly recommended that the surgeon perfect the technique for implantation of the TMJ Fossa-Eminence Prosthesis System™ and TMJ Condylar Prosthesis System™ through attendance at surgical demonstration courses, use of an instructional video, and manipulation of replica models. Instructional videos and literature are available from TMJ Implants, Inc. TMJ Implants, Inc. can provide names of individuals, independent from the company, with extensive experience for consultation prior to surgery.

- **TEST FOR ANY SUSPECTED SENSITIVITY TO MATERIALS**

Patients with suspected sensitivity to metals, such as Nickel, should undergo appropriate testing for sensitivity to Co-Cr-Mo alloy. Upon request, TMJ Implants, Inc. will supply a sample of this alloy and/or the chemical composition for pre-operative allergy testing. The device should not be used in patients who test positive for Co-Cr-Mo alloy sensitivity.

- **IF LONGER SCREWS ARE NECESSARY**

Occasionally, longer screws will be necessary to engage bone. It is important that the surgeon exercise great care to prevent injury to deeper vital structures. Care must be exercised not to penetrate or impinge any auditory structure, middle cranial fossa, or any neuro/vascular structures.

- **IF EXCISING BONE**

When performing an excision of bone in the area of the normal glenoid fossa and condyle, especially in cases of bony ankylosis, the surgeon must exercise great care to avoid penetration into the middle cranial fossa, the auditory canal, or other vital structures.

5. PRECAUTIONS

PRIOR TO SURGERY

- Special attention should be paid to patient selection. Careful evaluation should be made of patients with disorders that might interfere with their ability to comply with the limitations and precautions necessary to achieve beneficial outcome from this implant.
- All TMJ Metal-on-Metal Total Joint Replacement prostheses, screws, drills, and the Condylar sizers are provided sterile. Inspect sealed sterile package before opening. If seal is broken, do not use. Do not resterilize.
- Prior to use, the Fossa-Eminence Sizer Kit and Instrument Kit, or System, containing screwdrivers and Fossa-Eminence holders must be sterilized as outlined in Section 9.
- The surface of the device must remain clean and free of debris prior to implantation.
- The implant must be handled only with talc-free gloves to avoid introduction of talc into the implantation site.
- The prostheses must be protected from scratching or bending prior to and during surgical implantation, as such damage may cause weakening or fatigue of the metal or fracture of the part.

DURING SURGERY

- The TMJ Prosthesis System must be secured only through the use of the drills and screws supplied by TMJ Implants, Inc. The screws and drills used with the TMJ Prosthesis System have been specifically selected by size to ensure correct fixation of each prosthesis when used as directed. Any use of substitute drill bits or screws not supplied by TMJ Implants, Inc. in the TMJ Prosthesis System may result in less than optimal long-term results and may adversely affect the performance of the prosthetic device.
- It is strongly recommended that at least four (4) Fossa-Eminence screws be used where practical to achieve firm fixation of the TMJ Fossa-Eminence prosthesis.
- It is strongly recommended that at least six (6) condylar screws for the Universal Condylar prosthesis and nine (9) condylar screws with the Christensen-Chase Condylar prosthesis be used, where practical, to achieve firm fixation of the TMJ Condylar prosthesis. Care must be taken to secure at least 3 screws in the topmost holes, where practical.
- It is recommended that the head of the TMJ Condylar prosthesis be placed into the prosthetic glenoid fossa portion of the TMJ Fossa-Eminence prosthesis with all interposing soft tissue removed. The Condylar prosthesis articulating surface should preferably be centered in the Fossa and should not contact the screws of the Fossa-Eminence prosthesis.

6. ADVERSE EVENTS

Adverse events observed in the clinical use of the TMJ Implants, Inc. Metal-on-Metal Total Joint Replacement System are listed below (in descending order of frequency of occurrence.)

- Postoperative pain, swelling, jaw muscle spasm (10 subjects)
- Facial nerve and muscle weakness or paralysis (9 Subjects)
- Dislocation of the joint (2 subjects)
- Infection (2 subjects)
- Degenerative joint changes and development of adhesions (2 subjects)
- Nausea and vomiting (2 subjects)
- Perioperative bleeding (1 subject)
- Seizures (1 subject)
- Malocclusion (1 subject)

In addition to the adverse events identified above, potential adverse events and complications associated with temporomandibular joint surgery and reconstruction may require further treatment and include but are not limited to:

- Hematoma formation
- Hemorrhage
- Dental malocclusion, jaw dysfunction, limited range of motion
- Foreign body or allergic reactions to the device materials
- Rejection of the device
- Wear, displacement of the device or implant loosening
- Fracture of the device
- Hearing problems
- Surgical damage to anatomical structures adjacent to the TMJ
- Patient discomfort
- Speech problems
- Facial deformity

7. CLINICAL DATA

Two studies were conducted that support the safety and effectiveness the TMJ Implants, Inc. Metal-on-Metal Total Joint Replacement System. The first is the TMJ Implants, Inc. Registry, the second a prospective clinical study.

The objectives of these studies were to demonstrate that the TMJ Implants, Inc. Metal-on-Metal Total Joint Replacement System [TMJ ArthroChrome condylar prosthesis used in conjunction with a TMJ Fossa-Eminence Prosthesis (total joint replacement)] significantly reduces TMJ pain and improves interincisal opening.

An additional objective of the Prospective Study included a review of the incidence of device related adverse events occurring during the study.

For each study, pain measurements were recorded using a 10cm Visual Analog Scale (VAS). The left side of the scale represented no pain while the right side represented the most severe pain imaginable. The patients were instructed to mark a vertical line on the scale to indicate their perceived level of pain. Interincisal opening was measured in millimeters using a Therabite™ Scale. The interincisal opening was measured at the point at which the patient cannot open his/her mouth any wider.

TARGET POPULATION

The target population for a TMJ metal-on-metal total joint replacement utilizing the TMJ Condylar prosthesis in conjunction with the TMJ Fossa-Eminence prosthesis are those patients suffering from deficiencies of the natural condyle in cases of serious adhesion, condylar destruction, ankylosis, avascular necrosis, intrinsic bone disease (neoplasia), congenital disease involving the TMJ, rheumatoid arthritis, osteoarthritis, foreign body giant cell reaction, previous failed implant surgery, or other pathology with resultant occlusal or functional deficiency.

TMJ Implants, Inc. Registry

The Registry is a collection of data on patients that receive a TMJ Implants, Inc. device. The primary purpose of the Registry is for device tracking as required by US Regulation, 21 CFR Part 821. Monitoring clinical data was a secondary function of the Registry in the past. Clinical data was gathered by asking the operating surgeon to voluntarily submit baseline assessments of pain, and interincisal opening and then was asked to submit clinical follow up reports at 6 months, 1 year, 18 months, 2 years, and yearly thereafter.

DEMOGRAPHICS

There were 425 total joint recipients representing 1309 devices included in the cross-section data set and 63 patients representing 204 devices in the cohort data set. The majority of the patients in each group received "stock" devices. There were 267 patients (63%) with stock devices in the cross-section group, of which 25 (6%) were of the Christensen/Chase model. There were 36 patients (57%) in the cohort group with stock devices of which 10 (16%) were of the Christensen/Chase model. The mean age of the cross-section group of patients was 42±12 years and 41±12 years for the cohort. There were 89% female in the cross-section and 91% female in the cohort. No difference between age and gender among the cross-section group and cohort was demonstrated, p>0.05.

REDUCTION IN PAIN AND DIET RESTRICTION: TOTAL JOINT RECONSTRUCTION

From the cross-section data set, there is a marked reduction in pain and diet restriction within the first month after surgery, as demonstrated by Figure 1. Patients appear to reach their greatest relief within 6 months of surgery and maintain that level of improvement or slightly above through 4 years implant duration. A similar trend of pain reduction is demonstrated with a cohort of 63 patients with complete data through 2 years implant duration, Table 2. Approximately one-third of the cohort has complete data at 3 years. These data are not included in the analysis. By applying a test of contrasts to the 2-year cohort, a

significant reduction in pain at every time period is demonstrated when compared to the pre-op value, p<0.0001. Additionally, an ANOVA F-Test was applied to the cohort. The overall test is significant at p<0.0001.

Figure 1: Total Joint, Reduction in Pain and Diet Restriction


Table 1: Total Joint, Pain and Diet Cross Section Data

Months	0	1	6	12	18	24	36	48
Mean	7.8	3.9	3.0	2.9	3.7	3.6	3.7	4.2
n	425	397	272	185	105	81	34	11

Table 2: Total Joint, Pain Cohort Data

Pain	0	1	6	12	24	36
8.0	4.4	3.2	3.1	3.4	3.3	
n	63	63	63	63	63	20

Diet	0	1	6	12	18	24	36	48
Mean	7.7	4.4	2.9	2.7	3.2	3	3.3	4.1
n	386	363	258	180	105	78	34	11

Improvement in Interincisal Opening: Total Joint Reconstruction

From the cross-section data set, there is a marked improvement in the interincisal opening at 6 months after surgery, as demonstrated by Figure 2. Patients appear to reach their greatest improvement at 6 months after surgery and maintain that level of improvement through 4 years implant duration. Due to the small sample size at 4 years, the slight decrease at this time period is thought not to be significant. A similar trend of improvement in opening is demonstrated with a cohort of 57 patients with complete data through 2 years implant duration, Table 4. Approximately one-third of the cohort has complete data at 3 years. By applying a Test of Contrasts to the 2-year cohort, a significant improvement in opening is demonstrated at every time period when compared to the pre-op value, p<0.0001. Additionally, an ANOVA F-Test was applied to the cohort. The overall test is significant at p<0.0001.

Figure 2: Total Joint, Improvement in Opening


Table 3: Total Joint, Opening Cross-Section Data

Opening	0	1	6	12	18	24	36	48
Mean	20.9	25.8	30.4	31.2	30.4	31.1	30.4	24.3
n	370	343	238	163	98	78	34	10

Table 4: Total Joint, Opening Cohort Data

Opening	0	1	6	12	24	36
21.2	25.1	31.4	32.8	32.4	32.5	
n	57	57	57	57	57	20

Prospective Study

Demographic data is available from 43 subjects receiving a metal-on-metal total joint replacement. The mean age of this population is 44 ±13 years with a range of 24 to 74 years. The majority, 40 (93%), are female and Caucasian, 37 (86%). There are two (2) African-Americans and two (2) Hispanics, and one (1) Native-American. The most frequently reported indications for total joint replacement are previously failed implant surgery, 31%; recurrent fibrous or bony ankylosis, 29%; degenerative joint disease, 13%; and trauma, 9%. The remaining 18% included resorptive joint pathology, inflammatory arthritis, internal derangement, and unreported.

REDUCTION IN PAIN AND DIET RESTRICTION: TOTAL JOINT RECONSTRUCTION

On average, for patients with total joint replacements, there has been a 64% reduction in pain and diet restriction by the 6-month post-op visit, Figure 3. Between 12 and 18 months, however, the mean pain and diet scores appear to increase generally above 3.0cm and then return to near the 6-month level by 24 months. A comparison of the historical data at 12 months from patients with pain scores >3.0cm and those with pain scores ≤3.0cm indicates that all 10 patients presenting with pain scores >3.0cm are all multiply operated patients and 7 of the 9 patients presenting with pain scores ≤3.0cm are multiply operated. However, of the former group (>3.0cm), besides being multiply operated, 6 of the 10 also had previously failed implant surgery, whereas of the latter group (≤3.0cm), only 1 of the 7 patients that are multiply operated had experienced failed implant surgery. All other confounding variables associated with the patients histories such as trauma, degenerative disease, ankylosis, osteochondritis dissecans or avascular necrosis, are of similar frequency among both groups. Therefore, it appears that those patients who have experienced previously failed implant surgery have a less stable post-op course than those patients who did not. A similar anomaly in the data at 18 months is also seen with the registry data for total joint

replacement, but not as marked, Figure 4. The results reported to date suggest that multiply operated patients may not experience as great a reduction in pain. The results are interim results and final conclusions cannot be drawn until all patients have completed at least 3 years and the data are collected and analyzed.

Figure 3: Reduction in Pain, Total Joint Replacement


Table 5: Mean values, Reduction in Pain, Total Joint


Pain	0	0.5	3	6	12	18	24	36
Months	0	0.5	3	6	12	18	24	36
Mean	8.2	6.0	3.0	2.9	4.3	5.6	3.5	3.4
N	42	39	29	22	19	8	10	2

Table 6: Mean values, Reduction in Diet, Total Joint

Diet	0	0.5	3	6	12	18	24	36
Months	0	0.5	3	6	12	18	24	36
Mean	8.6	6.8	3.6	2.7	3.4	4.3	2.9	1
N	42	39	29	22	19	8	10	2

Note: There were only 42 pre-op pain and opening measurements from 43 patients. One patient failed to record a pre-op pain or opening measure; however, pain and opening were subsequently reported post-op through 18 months follow up for this patient.

Figure 4: Pain Reduction, Total Joint, Clinical Study and Registry Data

IMPROVEMENT IN INTERINCISAL OPENING: TOTAL JOINT RECONSTRUCTION

A patient's interincisal opening is measured in millimeters using either a Therabite scale or E-Z Flex system. Measurements are taken at the pre-op visit, within 10 days post-op, then at 3, 6, 12, 18, 24, and 36 months after implant. On average, by the 3-month post-op visit there is a significant improvement in opening, Figure 5. These trends are similar to the results presented with the TMJ Registry data, from those patients implanted with a total joint prosthesis, Figure 6. This trend of improvement continues through 36 months post-implant, Table 7.

Figure 5: Improvement in Opening, Total Joint replacement



Table 7: Opening Mean values, Total Joint replacement

Opening	0	0.5	3	6	12	18	24	36
Months	0	0.5	3	6	12	18	24	36
Mean	21.9	22.1	30.3	31.8	31.8	27.4	33.4	40.5
N	37	40	30	22	19	9	12	2

Figure 6: Improvement in Opening, Total Joint, Clinical Study and Registry Data

8. INFORMATION FOR USE

The TMJ Metal-on-Metal Total Joint Replacement System™ is comprised of the TMJ Fossa-Eminence Prosthesis™ and the TMJ Condylar Prosthesis™.

This Instructions for Use is intended to give you some answers for the use of TMJ Implants, Inc.'s TMJ Metal-on-Metal Total Joint Replacement System™. However, this is not intended to be an exhaustive, or comprehensive treatise on this subject of alloplastic joint reconstruction.

PRIOR TO SURGERY

Warnings

Patients with suspected sensitivity to metals, such as Nickel, should undergo appropriate testing for sensitivity to Co-Cr-Mo alloy. Upon request, TMJ Implants, Inc. will supply a sample of this alloy and/or the chemical composition for pre-operative allergy testing. The device should not be used in patients who test positive for Co-Cr-Mo alloy sensitivity.

There are instances where this technique is not recommended due to prior surgical procedures and the need to place prostheses in less than optimal angles and positions, or in cases where systemic medical disease would contraindicate this implant procedure in the view of the operating surgeon. The operating surgeon must make this evaluation. It is the surgeon's responsibility to determine the need for patient specific implants given anatomical considerations or unusual masticatory forces in a given patient.

Patients undergoing local or general anesthesia, prolonged dental therapy, extraction of teeth, or those patients using mechanical devices which create abnormal forces within the joint need to be alerted to possible injury to the joint or prosthesis due to those unusual forces.

The TMJ Metal-on-Metal Total Joint Replacement System™ is intended for total joint replacement. The TMJ Condylar prosthesis must be implanted in conjunction with a TMJ Fossa-Eminence prosthesis manufactured by TMJ Implants, Inc.

All TMJ Metal-on-Metal Total Joint Replacement prostheses, screws, drills, and the Condylar sizers are provided sterile.

Prior to use, the Fossa sizer, and the screwdrivers and holders for both prostheses must be sterilized in their respective containers in accordance with hospital standards for steam sterilization (see Section 9).

Caution - The surface of the device must remain clean and free of debris prior to implantation.

Caution - The implant must be handled only with talc-free gloves to avoid introduction of talc into the implantation site.

Caution - The prostheses must be protected from scratching or bending prior to and during surgical implantation, as such damage in certain circumstances may cause weakening or fatigue of the metal or fracture of the part.

A detailed Patient-Specific Manual is available which provides instructions for CT scanning patients, and describes how to prepare the Anatomical Bone Model prior to implant design.

It is the responsibility of the surgeon to become familiar with the surgical techniques for implantation of these devices through attendance at surgical demonstration courses, use of instructional video, consultation with experienced associates and manipulation of replica models.

DURING SURGERY

The location of each component included in the TMJ Metal-on-Metal Total Joint Replacement Prosthesis System™ is shown on a diagram on the inside lid of the respective package.

The TMJ Fossa-Eminence prosthesis used in conjunction with the TMJ Condylar prosthesis should be placed and secured first.

PLACING THE TMJ FOSSA-EMINENCE PROSTHESIS COMPONENT OF THE TMJ METAL-ON-METAL TOTAL JOINT REPLACEMENT SYSTEM™

The normal preauricular or endaural incision and approach to the joint is accomplished. Exposure of the entire zygomatic process of temporal bone lateral to the joint is necessary to facilitate placement of the TMJ Fossa-Eminence prosthesis.

When the joint is fully exposed try the sizer for fit. Take your time at this point. Find the sizer that fits the bone most accurately with at least 3-point contact and allows the condyle to function smoothly, without dislocation of the joint and provides suitable stability. At this stage check the occlusion very carefully. Ensure the occlusion remains as seen pre-operatively or as desired post-operatively. If not, determine why.

After selecting the proper sizer, check the laser-etched number on the sizer and have the nurse or anesthesiologist record it for future reference.

Noting the etched number from the correct-fitting sizer, select the same numbered TMJ Fossa-Eminence prosthesis, which has been packaged sterile. Try it for accuracy of fit, proper occlusion, and mobility of the condyle.

It is strongly recommended that at least four (4) Fossa-Eminence screws be used where practical to achieve firm fixation of the TMJ Fossa-Eminence prosthesis. Be sure to use the screws provided to insure compatibility of the metals. Caution should be used so as not to force the screw in place with too much pressure as the screw head could fracture. Always drill the hole slightly deeper than the length of the screw. When the implant has been secured in place, check again for proper jaw function and proper occlusion. Be sure to use the drill bits provided for preparing the screw holes. Be diligent in this surgery to avoid injury to important adjacent structures i.e. middle cranial fossa, ear structures, facial nerve, and middle meningeal artery.

The TMJ Fossa-Eminence prosthesis must be secured only through the use of the drills and screws supplied by TMJ Implants, Inc. The screws and drills used with the TMJ Fossa-Eminence prosthesis have been specifically selected by size to ensure correct fixation of each prosthesis when used as directed. Any use of substitute drill bits or screws not supplied by TMJ Implants, Inc. in the TMJ Fossa-Eminence Prosthesis System™ may result in less than optimal long-term results and may adversely affect the performance of the prosthetic device.

PLACING THE TMJ CONDYLAR PROSTHESIS COMPONENT OF THE TMJ METAL-ON-METAL TOTAL JOINT REPLACEMENT SYSTEM™

A CT Scanning Protocol is available directly from TMJ Implants, Incorporated, upon request.

It is the responsibility of the surgeon to become familiar with the surgical techniques for implantation of these devices through attendance at surgical demonstration courses, use of instructional video, consultation with experienced associates and manipulation of replica models.

Once the wound has been opened to expose the remaining natural condyle and the occlusion have been fixed, then remove enough condyle height to allow a condylar sizer with its 13mm head to be placed in position. The condylar sizer will allow you to determine which length of the TMJ Condylar prosthesis will most accurately fit the patient's mandible.

When the correct condylar length has been determined by use of the condylar sizer, you are now ready to place the actual TMJ Condylar prosthesis.

It is recommended that the head of the TMJ Condylar prosthesis be placed into the prosthetic glenoid fossa portion of the TMJ Fossa-Eminence prosthesis with all interposing soft tissue removed. The Condylar prosthesis articulating surface should preferably be centered in the Fossa and should not contact the screws of the Fossa-Eminence prosthesis.

The proper placement of the condylar head into the TMJ Fossa-Eminence prosthesis, assures that the head does not contact any screw heads during function. It is important to fix the TMJ Condylar prosthesis to the ramus of the mandible with as many screws as possible. It is strongly recommended that at least six (6) condylar screws for the Universal Condylar prosthesis and nine (9) condylar screws with the Christensen-Chase Condylar prosthesis be used, where practical, to achieve firm fixation of the TMJ Condylar prosthesis. Care must be taken to secure at least 3 screws in the topmost holes, where practical.

The Christensen/Chase TMJ Condylar Prosthesis System is available in either a left or right-sided set. The head of each of these prostheses is identical in size and shape. As with the Universal System, there are three lengths available: 45mm, 50mm, and 55mm. The Christensen/Chase prostheses differ from the Universal prostheses in that the Christensen/Chase devices have an angled extension on the distal portion of the flange, allowing the physician to more closely follow the patient's natural mandibular structure, and to provide anchoring options in the absence of bone.

The flange direction of the Condylar prosthesis is generally ideal when it parallels the posterior margin of the mandibular ramus.

The TMJ Condylar prosthesis must be secured only through the use of the drills and screws supplied by TMJ Implants, Inc. The screws and drills used with the TMJ Condylar prosthesis have been specifically selected by size to ensure correct fixation of each prosthesis when used as directed. Any use of substitute drill bits or screws not supplied by TMJ Implants, Inc. in the TMJ Condylar Prosthesis System may result in less than optimal long-term results and may adversely affect the performance of the prosthetic device.

If both sides of the mandible area are to be reconstructed, the contralateral side is operated in a similar fashion.

The mandible should now be mobilized and forced open, observing the maximum range of opening. Should mobility be limited, the surgeon should assess the possible need for coronoidectomies and/or relief of adhesions to the ramus. Observe the movement of the condylar head(s) to be sure that dislocation does not occur. The fixation appliances may now be removed.

ADDITIONAL CONSIDERATIONS

The placement of a fat graft around the implants at surgery may reduce the occurrence of subsequent adhesion or even ankylosis. For those patients susceptible to heterotopic bone formation, appropriate fat grafts and radiation therapy should be considered.

Occasionally, longer screws will be necessary to engage bone. It is important that the surgeon exercise great care to prevent injury to deeper vital structures. Care must be exercised not to penetrate or impinge any auditory structure, middle cranial fossa, or any neuro/vascular structures.

When performing an excision of bone in the area of the normal glenoid fossa and condyle, especially in cases of bony ankylosis, the surgeon must exercise great care to avoid penetration into the middle cranial fossa, the auditory canal, or other vital structures.

POST SURGERY

Accepted surgical practice should be followed in post-operative care.

After the total joint replacement is completed, all instruments must be thoroughly cleaned, decontaminated and sterilized in accordance with the following procedures outlined below.

9. CLEANING AND STERILIZATION

DESCRIPTION	Provided	Use
TMJ Fossa-Eminence Prosthesis System™	Sterile	Single
TMJ Fossa-Eminence Trial Sizing System	Non-Sterile	Reusable
TMJ Condylar Prosthesis System™	Sterile	Single
TMJ Condylar Trial Sizers	Sterile	Single
TMJ Christensen/Chase Condylar Prosthesis System	Sterile	Single
TMJ Christensen/Chase Condylar Prostheses Sizers	Sterile	Single
TMJ Patient-Specific™ Fossa-Eminence Prostheses	Sterile	Single
TMJ Patient-Specific™ Condylar Prostheses	Sterile	Single
TMJ Screw and Drill Kit	Sterile	Single
TMJ Instrument Kit	Non-Sterile	Reusable
TMJ Fossa-Eminence Trial Sizing System & Instrument System	Non-Sterile	Reusable

STERILE COMPONENTS PROVIDED

Contents must be stored at 10 to 32 degrees C (50-90° F) at 20%-80% relative humidity.

The TMJ Condylar Prostheses, condylar sizers, screws, and drills, and the Fossa-Eminence prostheses, screws and drills are all packaged sterile in double Tyvek/film or paper/film peel pouches utilizing gamma irradiation or electron beam (e-beam) irradiation sterilization at levels of 25 - 40 kGy with a sterility assurance level (SAL) of 10⁻⁶. Care must be taken to assure packaging remains undamaged to ensure sterility.

CLEANING INSTRUCTIONS FOR REUSABLE COMPONENTS

For your safety, be familiar with the procedures for handling contaminated materials at your facility prior to utilizing these instructions.

Clean instruments before sterilization in the autoclave trays provided by TMJ Implants, Inc. as soon as possible **after** use. Avoid allowing soiled instruments to dry. Immerse into or use towels dampened with deionized or distilled water to keep soiled instruments moist prior to cleaning.

Manually wash the sizer templates and instrumentation with mild enzymatic detergent following the manufacture's instructions for use. Follow the detergent manufacturer's recommendations for use dilution. Use enzymatic detergents (e.g. Enzol™) and warm/hot water prepared as recommended by the manufacturer. PH neutral detergents are recommended. If acidic or alkaline solutions are used, follow the manufacturer's recommendations for neutralizing the pH by rinsing with water or other neutralizing solution. Highly alkaline cleaners (pH ≥ 12) are not recommended. Avoid exposure to acidic or alkaline solutions and solutions containing chlorides, bromides, or iodines.

Allow the devices to soak in a mild enzymatic detergent for a minimum of one minute. Use a soft bristle brush to manually clean the devices while immersed in the cleaning detergent, paying particular attention to crevices and other hard-to-clean areas. Clean the devices until all adherent visible soil is removed.

After washing, thoroughly rinse instruments for one minute under lukewarm clean, deionized or distilled water.

Inspect for cleanliness, especially in crevices or recesses. Check instruments thoroughly for damage, (i.e. chips, cracks, corrosion, surface wear, etc.) especially with instruments with moving parts or interfits such as a quick-connect mechanism. Do not use instruments that have been damaged. Damaged instruments should be replaced before further use.

Dry completely with a clean, soft cloth before sterilization.

STERILIZATION INSTRUCTIONS FOR REUSABLE COMPONENTS

Reusable instruments, e.g., the Fossa sizer, and the screwdrivers and holders, for both prostheses must be sterilized prior to surgical use and after cleaning post-surgically in their respective containers provided by TMJ Implants, Inc.

The following sterilization cycle has been shown to sterilize product to a sterility assurance level (SAL) of 10⁻⁶. Other similar cycles may be used but have not been evaluated. It is the responsibility of the customer to demonstrate the appropriateness of the sterilization cycle used should it vary from the following:

Sterilizer Type:	Prevacuum Steam Sterilization
Preconditioning Pulses:	3
Minimum Temperature:	132° C
Full Cycle Time:	4 minutes
Minimum Dry Time:	20 minutes
Maximum Load Configuration:	5 trays
Minimum Load Configuration:	1 tray
Sample Configuration:	Wrapped

It is recommended not to exceed 5-stacked trays per sterilization run.

LIMITED WARRANTY

TMJ Implants, Inc. warrants that this product meets the manufacturer's specifications and is free from manufacturing defects at the time of delivery. This warranty specifically excludes defects resulting from misuse, abuse, or improper handling of the product subsequent to receipt by the purchaser.

**Should you have any questions, please feel free to call us at
(303) 277-1338 or (800) 825-4865.**

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